

Email :
sales@schooleducationalinstruments.com
Phone: +91-0171-2601773

Product Name :
Geared Teacher Clock for School Maths Lab

Product Code :
MathsLab0120

Description :

Geared Teacher Clock This Geared clock Accurate to the hours with the movement of minute Arm. This geared demonstration clock is a great manipulative of time-telling skills for your students. This clock is provided with bright color clock arms, featuring easy-to-read. Hidden gears maintain correct hour and minute relationships as you manipulate movable hands. Clock is made of durable plastic.

NAUGRA[®]

School Educational Instruments

Website: [www.https://www.schooleducationalinstruments.com/](https://www.schooleducationalinstruments.com/), **Email:** sales@schooleducationalinstruments.com

Address: Ambala Cantt, Haryana, India **Phone:** 91-0171-2643080