

Email :  
sales@schooleducationalinstruments.com  
Phone: +91-0171-2601773

**Product Name :**  
Geared Student Clock for School Maths Lab

**Product Code :**  
MathsLab0123

**Description :**

Geared Student Clock This Student clock Accurate to the hours with the movement of minute Arm. This geared student clock is a great manipulative to understand time-learning skills for individual student. This clock is provided with bright color clock arms, featuring easy-to-read. Hidden gears maintain correct hour and minute relationships as you move pully on the top of the clock. This Clock is made of durable plastic.

**NAUGRA<sup>®</sup>**

**School Educational Instruments**

**Website:** [www.https://www.schooleducationalinstruments.com/](https://www.schooleducationalinstruments.com/), **Email:** sales@schooleducationalinstruments.com

**Address:** Ambala Cantt, Haryana, India **Phone:** 91-0171-2643080