

Email :
sales@schooleducationalinstruments.com
Phone: +91-0171-2601773

Product Name :
Cement Samplers

Product Code :
SCHOOL-CEM740001



Description :

Cement Samplers

Technical Specification :

Made of brass, it has an outside dia. of 32 mm x 700 mm long. Used to sample cement from packages. Weight approx.: 2 kg Bulk cement sampler It consists of two brass concentric tubes with slots. The inner tube rotates to close the slots and take the sample. Used to sample cement in bulk storage or bulk shipments. The capacity of the internal tube is approx. 3 litres. Weight approx.: 1 kg

NAUGRA[®]

School Educational Instruments

Website: [www.https://www.schooleducationalinstruments.com/](https://www.schooleducationalinstruments.com/), **Email:** sales@schooleducationalinstruments.com

Address: Ambala Cantt, Haryana, India **Phone:** 91-0171-2643080