

Email :
sales@schooleducationalinstruments.com
Phone: +91-0171-2601773

Product Name :
Le Chatelier Moulds

Product Code :
SCHOOL-CEM740003



Description :

Le Chatelier Moulds

Technical Specification :

The mould consists of a spring tensioned split cylinder 30 mm internal diameter, 30 mm high with two indicator stems which measure 165 mm from the points of the centre line of the cylinder and O ring. Used for determining the expansion of cement. Two or three moulds are required for each test. To perform the test, a water bath is also required. See accessories. We propose three packages : Le Chatelier Mould Single unit. Individually checked, supplied complete with certificate of conformity and identified by a serial number. Weight approx.: 30 g Le Chatelier Mould. Pack of 6 units. Weight approx.: 150 g Le Chatelier Soundness Kit The kit includes all accessories to perform the test and to verify the conformity of the moulds. Consisting of : Le Chatelier Moulds Glass plates, 50 mm sq.

NAUGRA[®]

School Educational Instruments

Website: [www.https://www.schooleducationalinstruments.com/](https://www.schooleducationalinstruments.com/), **Email:** sales@schooleducationalinstruments.com

Address: Ambala Cantt, Haryana, India **Phone:** 91-0171-2643080