

Email :
sales@schooleducationalinstruments.com
Phone: +91-0171-2601773

Product Name :
Microprocessor PH Meter

Product Code :
SCHOOL-EVN750002



Description :

Microprocessor PH Meter

Technical Specification :

Display: 16 X 2 Alpha Numeric LCD Display. Accuracy: ± 0.01 pH, $\pm .1$ mV, Temp. $\pm 1^\circ\text{C}$. (Auto/Manual). Temp. Compensation: 0 to 100°C . Range: pH : 0 to 14.00, mV : 0 to ± 1999.9 . Resolution: 0.01 pH , 0.1 mV, Temp. 0.1°C . Calibration: Auto / Manual. Storage: Up to 100 Samples. Power: 230V $\pm 10\%$, 50Hz AC. Auto Buffer Recognition: Auto 4.00, 7.00, 9.20 pH.

NAUGRA[®]

School Educational Instruments

Website: [www.https://www.schooleducationalinstruments.com/](https://www.schooleducationalinstruments.com/), **Email:** sales@schooleducationalinstruments.com

Address: Ambala Cantt, Haryana, India **Phone:** 91-0171-2643080