

Email :
sales@schooleducationalinstruments.com
Phone: +91-0171-2601773

Product Name :
Basic Refrigeration and Air Conditioning Training
System

Product Code :
SCHOOL-REF870002



Description :

Basic Refrigeration and Air Conditioning Training System

Technical Specification :

The Basic Refrigeration And Air Conditioning Training System main components are compressor, condenser and receiver plus electrical and communications systems. The models are plugged onto the base unit and connected hydraulically with refrigerant hoses and electrically with cables. Self-sealing couplings reduce the refrigerant loss to a minimum. All components are arranged well visible to allow their operation to be monitored. The modern and powerful software is an integral part of the training system in the form of hardware/software integration (HSI). It enables the comfortable execution and analysis of the experiments. Temperatures and pressures in the system are recorded by sensors and displayed dynamically in the software for system operation and data acquisition. The experimental unit is connected to the PC via a USB interface. FEATURES: Fundamentals of the refrigeration cycle Fundamentals of air conditioning Components in a refrigeration system/air conditioning system System operation Fault finding Base unit for the setup of basic experiments in refrigeration and air conditioning technology SPECIFICATIONS: Condensing unit Refrigeration capacity: 340W at 0/32°C Receiver: 0,7L Refrigerant currently available in the market Filling volume: 700g CO2-equivalent: 0,4t Measuring ranges : Temperature: 1x -50...50°C, 3x 0...100°C Pressure : 1x intake side: -1...9bar 2x delivery side: -1...15bar Required for operation : 230V, 50Hz, 1 phase 230V, 60Hz, 1 phase

NAUGRA[®]

School Educational Instruments

Website: [www.https://www.schooeducationalinstruments.com/](https://www.schooeducationalinstruments.com/), **Email:** sales@schooeducationalinstruments.com

Address: Ambala Cantt, Haryana, India **Phone:** 91-0171-2643080