

Email :  
sales@schooleducationalinstruments.com  
Phone: +91-0171-2601773

**Product Name :**  
Say Bolt Viscometer Apparatus

**Product Code :**  
SCHOOL-FLE910001



**Description :**

Say Bolt Viscometer Apparatus

**Technical Specification :**

The Saybolt Viscometer is used to determine empirical measurement of Saybolt Viscosity of petroleum products at specified temperatures. The time it takes for a 60mL sample to flow through a tube with a specified diameter (furol or universal) is measured and reported in seconds. The viscometer can be used for temperatures between 21 and 99 °C (70 to 210 °F). The viscometer includes water-oil bath, stirrer, cooling coil, electric heater with digital thermo regulator. Viscosity Thermometer set consists of 6 thermometers with the temperature ranges; 19 to 27°C, 34 to 42°C, 49 to 57°C, 57 to 65°C, 79 to 87°C (250 mm length) and 95 to 103°C where each thermometer has 0.1°C subdivisions. Filter funnel, withdrawal tube and thermometer set should be ordered separately. The Saybolt Two-Tube Digital Viscometer is supplied complete with: Universal Orifice, 2 pcs Furol Orific, , 2 pcs Heat Transfer Oil, 5 lt Key Saybolt viscosity flask, glass, 60 ml, 2 pcs. Thermometer Support, 2 pcs Stopper, 2 pcs Technical Specifications : Dimensions: 450x300x550 mm Weight (approx.): 10 kg Power: 750 W

**Address:** Ambala Cantt, Haryana, India **Phone:** 91-0171-2643080

**NAUGRA<sup>®</sup>**