

Email :
sales@schooleducationalinstruments.com
Phone: +91-0171-2601773

Product Name :
Freezing Point Apparatus

Product Code :
SCHOOL-FLE910018



Description :

Freezing Point Apparatus

Technical Specification :

Double tube Ø 30 x 240 mm fitted with cap

Dewar jar Ø 75 x 280 mm

Geared motor for stirring at 80rpm with wire stirrer

Digital thermometer

Mounted on a plate structure painted with anti-acid epoxidic products

Freezing Point Apparatus

ACCESSORIES:

ASTM-114C Thermometer ASTM 114C IP 14C

102-376 Type A Nitrogen collar made in glass

SPARE PARTS:

102-371 Double tube 30 x 240mm

102-372 Wire stirrer

102-373 Dewar jar

102-374 Cork series

102-375 Type B Modified moisture proof collar made in glass

Electrical:

specify 120 or 230 Volts, 50/60 Hz

NAUGRA[®]

School Educational Instruments

Website: [www.https://www.schooleducationalinstruments.com/](https://www.schooleducationalinstruments.com/), **Email:** sales@schooleducationalinstruments.com

Address: Ambala Cantt, Haryana, India **Phone:** 91-0171-2643080