

Email :
sales@schooleducationalinstruments.com
Phone: +91-0171-2601773

Product Name :
Aneroid Barometer

Product Code :
SCHOOL-MES880039



Description :

Aneroid Barometer

Technical Specification :

FEATURES:

Bimetal element provides temperature compensation for aneroid capsule.

Choose from mbar/mm Hg or mbar/inches Hg models.

Includes an attractive brass case and include a bezel for easy bulk-head mounting.

Needs no batteries or external power.

Use these economical instruments for long-term temperature, relative humidity and barometric pressure monitoring.

SPECIFICATION:

Triple Scale Aneroid Barometer is made of solid brass & stainless steel construction, which keeps it durable for years of classroom use.

Measuring 6 in diameter it is wall mountable, calibrate able and is also RoHS compliant.

The range of the meter is 27.75 – 31.25 in Hg, 70.5 – 79.4 cm Hg, 940 – 1080 m bar and the accuracy of the meter is ± 0.1 in Hg, ± 0.3 cm Hg, ± 3.5 m bar.

NAUGRA[®]

School Educational Instruments

Website: [www.https://www.schooleducationalinstruments.com/](https://www.schooleducationalinstruments.com/), **Email:** sales@schooleducationalinstruments.com

Address: Ambala Cantt, Haryana, India **Phone:** 91-0171-2643080