

Email :
sales@schooleducationalinstruments.com
Phone: +91-0171-2601773

Product Name :
Digital Thermometer

Product Code :
SCHOOL-MES880041



Description :

Digital Thermometer

Technical Specification :

FEATURES:

High Resolution 0.01°C

Better Than 40 milliKelvin (0.04°C) System Accuracy

Digital Matching of Calibrated Sensors

Very High Accuracy and Stability

Circuitry Self-calibrates for Total Stability

3 or 4 Wire Sensors with Automatic Recognition

Readout Directly in °C, °F, Kelvin, and Ohms

Bench Mounting (Or Panel Option)

Rechargeable Battery/ ac Adapter

SPECIFICATION:

High-accuracy digital thermometer

Metal, bench-top; optional panel mounting kit

Precision Pt100 probe, temperature range -50 to 250°C (-58 to 482°F)

NAUGRA[®]

School Educational Instruments

Website: [www.https://www.schooeducationalinstruments.com/](https://www.schooeducationalinstruments.com/), **Email:** sales@schooeducationalinstruments.com

Address: Ambala Cantt, Haryana, India **Phone:** 91-0171-2643080