

Email :
sales@schooleducationalinstruments.com
Phone: +91-0171-2601773

Product Name :
Dial Gauge with Magnetic Stand

Product Code :
SCHOOL-MES880043



Description :

Dial Gauge with Magnetic Stand

Technical Specification :

FEATURES:

0.001" X 1" travel dial indicator.

Dial indicator crystal is hard coated to resist scratches.

Magnetic Stand accepts all AGD-type dial indicators, test indicators and probes having .375", 6mm or 8mm diameter stems and/or lug back with 1/4" diameter mounting hole.

Magnetic Stand has 132 Lbf of magnetic force to help prevent damage to indicators and probes due to accidental base fall-over.

The On/Off switch offers instant mounting and dismounting without any adverse effect to the indicators or work surface.

SPECIFICATION:

Magnetic base engages with just the turn of a switch and allows pinpoint adjustment.

The dial indicator features 0-1-Inch travel and has a resolution of 0.001-Inch.

Set includes a molded case for protection and convenience.

NAUGRA[®]

School Educational Instruments

Website: [www.https://www.schooeducationalinstruments.com/](https://www.schooeducationalinstruments.com/), **Email:** sales@schooeducationalinstruments.com

Address: Ambala Cantt, Haryana, India **Phone:** 91-0171-2643080