


Email :
sales@schooleducationalinstruments.com
Phone: +91-0171-2601773

Product Name :
Gas Charging Trainer

Product Code :
SCHOOL-REF870011



Description :

Gas Charging Trainer

Technical Specification :

FEATURES:

Evacuating a refrigeration system

Charging a refrigeration system

SPECIFICATION:

Portable bench top unit

Vacuum pump and filling cylinder

Manometer for intake pressure, delivery pressure and cylinder pressure; vacuum gauge

Manometer for intake and delivery pressure with temperature scale for refrigerant R134a

Connections for intake and delivery pressure

Filling cylinder: 1000g R134a

Measuring ranges:

Pressures intake: -1...10 bar

Delivery: -1...30 bar

Vacuum: 0...1000 mbar

Cylinder: -1...25 ba

Required for operation :

230V, 50Hz, 1 phase

230V, 60Hz, 1 phase

NAUGRA[®]

School Educational Instruments

Website: [www.https://www.schooleducationalinstruments.com/](https://www.schooleducationalinstruments.com/), **Email:** sales@schooleducationalinstruments.com

Address: Ambala Cantt, Haryana, India **Phone:** 91-0171-2643080