

Email :
sales@schooleducationalinstruments.com
Phone: +91-0171-2601773

Product Name :
Hydraulic Ram Test Bed

Product Code :
SCHOOL-FUD810018



Description :

Hydraulic Ram Test Bed

Technical Specification :

Pump body manufactured from clear acrylic with stainless steel pulse and non-return valves

Adjustable acrylic header tank with inlet and outlet hoses

Outlet hose with variable head arrangement

Demonstration of the water hammer effect to produce a pumping action

Quick-release fitting for easy connection to Hydraulics Bench

Supplied with weights to load pulse valve

Educational software available as an option

Supply head: 300-700mm variable

Delivery head: 750-1500mm variable

A supply pipe fitted with an inner and outer valve,

An air vessel to reduce hydraulic shock, and a delivery section.

NAUGRA[®]

School Educational Instruments

Website: [www.https://www.schooleducationalinstruments.com/](https://www.schooleducationalinstruments.com/), **Email:** sales@schooleducationalinstruments.com

Address: Ambala Cantt, Haryana, India **Phone:** 91-0171-2643080