

Email :
sales@schooleducationalinstruments.com
Phone: +91-0171-2601773

Product Name :
Computerized Reciprocating Pump Test Apparatus

Product Code :
SCHOOL-TUB820006



Description :

Computerized Reciprocating Pump Test Apparatus

Technical Specification :

FEATURES:

Principle of operation of a piston pump

Determination of efficiencies

Recording of pump characteristics

Pressure curves of delivery pressure and cylinder pressure

Influence of pulsation damping

SPECIFICATION:

Gear transmission ratio: 1=7,5

Piston pump :

Speed: 30...180min⁻¹

Max. Flow rate: 135L/h

Max. Head: 40m

Drive motor :

Power: 180W

Measuring ranges :

Pressure (cylinder): 0...5bar

Pressure (outlet): 0...5bar

Crank angle: 0...360°

Flow rate: 0,2...6L/min

230V, 50Hz, 1 phase

230V, 60Hz, 1 phase; 120V, 60Hz, 1 phase

NAUGRA[®]

School Educational Instruments

Website: [www.https://www.schooleducationalinstruments.com/](https://www.schooleducationalinstruments.com/), **Email:** sales@schooleducationalinstruments.com

Address: Ambala Cantt, Haryana, India **Phone:** 91-0171-2643080