

Email :
sales@schooleducationalinstruments.com
Phone: +91-0171-2601773

Product Name :
Separating and Throttling Calorimeter with Mini Boiler
Test Bench

Product Code :
SCHOOL-THM830023



Description :

Separating and Throttling Calorimeter with Mini Boiler Test Bench

Technical Specification :

FEATURES:

Determining the vapour content using

A separating calorimeter with cyclone water separator

A throttling calorimeter with vapour depressurization

Using an h-s diagram

SPECIFICATIONS:

Supplied vapour :

Temperature: 240°C

Pressure: 10bar

Safety valve: 10bar

Measuring ranges :

Temperature: 0...400°C

Pressure (outlet): -150...100mbar

Pressure (inlet): 0...16bar

230V, 50Hz, 1 phase

230V, 60Hz, 1 phase

120V, 60Hz, 1 phase

NAUGRA[®]

School Educational Instruments

Website: [www.https://www.schooleducationalinstruments.com/](https://www.schooleducationalinstruments.com/), **Email:** sales@schooleducationalinstruments.com

Address: Ambala Cantt, Haryana, India **Phone:** 91-0171-2643080