

Email :  
sales@schooleducationalinstruments.com  
Phone: +91-0171-2601773

**Product Name :**  
Power Door Lock System Trainer

**Product Code :**  
SCHOOL-ATU840021



**Description :**

Power Door Lock System Trainer

**Technical Specification :**

**FEATURES:**

OEM components (door lock actuators, switches, with OEM matching wire colors).

A lockable system fault box allowing instructors to insert from one to four faults for student diagnostic learning.

Electrical test points are included on the trainer for circuit diagnosis.

Students can wire this circuit using supplied leads and included OEM wiring diagrams.

Trainer requires an external 12V DC power supply to operate.

User manual included.

**SPECIFICATIONS:**

Power door lock actuators (2ea).

Door lock relay.

Fault box with 4 faults.

Test points for diagnosis.

Set of 21-jumper leads, 24" (60cm) (6 red, 1 black, 14 blue).

**NAUGRA<sup>®</sup>**

## **School Educational Instruments**

**Website:** [www.https://www.schooeducationalinstruments.com/](https://www.schooeducationalinstruments.com/), **Email:** [sales@schooeducationalinstruments.com](mailto:sales@schooeducationalinstruments.com)

**Address:** Ambala Cantt, Haryana, India **Phone:** 91-0171-2643080