

Email :
sales@schooleducationalinstruments.com
Phone: +91-0171-2601773

Product Name :
School Physics Light and Optical Lab Equipment
Catalog

Product Code :
SCHOOL-LGT230000



Description :

School Physics Light and Optical Lab Equipment Catalog

Technical Specification :

SCHOOL-LGT230001	School Semiconductor Laser Physics Lab
SCHOOL-LGT230002	School Magnifying Glasses With Led light
SCHOOL-LGT230003	School Spectrum Tube Power Supply Physics Lab
SCHOOL-LGT230004	School Concave Mirror With Stand
SCHOOL-LGT230005	School Concave Lens
SCHOOL-LGT230006	School Triangular Prism
SCHOOL-LGT230007	School Three Folding Magnifier
SCHOOL-LGT230008	School Handheld Reading Magnifier
SCHOOL-LGT230009	School Reading Magnifier
SCHOOL-LGT230010	School Pocket Magnifier

SCHOOL-LGT230011	School Handheld Magnifier
SCHOOL-LGT230012	School Optical Lens Set
SCHOOL-LGT230013	School Optical Experiment Set
SCHOOL-LGT230014	School Optical Tool Kit
SCHOOL-LGT230015	School Optical Bench
SCHOOL-LGT230016	School Optical Bench Set
SCHOOL-LGT230017	School Labs Optical Kit
SCHOOL-LGT230018	School Laser Ray Box
SCHOOL-LGT230019	School Physical Optical Kit
SCHOOL-LGT230020	School Tertiary Color Geometrical Optics Kit
SCHOOL-LGT230021	School Handheld Periscope
SCHOOL-LGT230022	School Handheld Stroboscope
SCHOOL-LGT230023	School Secchi Disk
SCHOOL-LGT230024	School Laser Refraction Tank
SCHOOL-LGT230025	School Physical Demonstration Polarizer
SCHOOL-LGT230026	School Mirage Magic Illusion
SCHOOL-LGT230027	School Plancks Constant Measuring Instrument
SCHOOL-LGT230028	School Kirchhoff Bunsen Spectroscope
SCHOOL-LGT230029	School Demonstrator of Polarization Application
SCHOOL-LGT230030	School He Ne leaser Instrument
SCHOOL-LGT230031	School Light Tricolor Experiment
SCHOOL-LGT230032	School Color Synthesis Experiment
SCHOOL-LGT230033	School Light Tricolor Demonstrator
SCHOOL-LGT230034	School Electric Newtons Color Disc
SCHOOL-LGT230035	School Hand Operated Newton Color Disk
SCHOOL-LGT230036	School Xenon Stroboscope

SCHOOL-LGT230037	School Spectral Lamp Power Supply with Lamp Holder
SCHOOL-LGT230038	School Optical Bench Kit
SCHOOL-LGT230039	School Sodium Spectral Lamp
SCHOOL-LGT230040	School Newtons Colour Disc
SCHOOL-LGT230041	School Polarising Filter Stress Demonstrator
SCHOOL-LGT230042	School Raybox Mirrored with Accessories
SCHOOL-LGT230043	School Direct Vision Spectroscope
SCHOOL-LGT230044	School Ray Box 12v
SCHOOL-LGT230045	School Prism Acrylic 22mm
SCHOOL-LGT230046	School Newtons Ring Apparatus
SCHOOL-LGT230047	School Red Laser Pen
SCHOOL-LGT230048	School Light Guide Demonstrator
SCHOOL-LGT230049	School Lens Holder Wooden

NAUGRA[®]

School Educational Instruments

Website: [www.https://www.schooleducationalinstruments.com/](https://www.schooleducationalinstruments.com/), **Email:** sales@schooleducationalinstruments.com

Address: Ambala Cantt, Haryana, India **Phone:** 91-0171-2643080