

Email :
sales@schooleducationalinstruments.com
Phone: +91-0171-2601773

Product Name :
20MHz Dual channel oscilloscope

Product Code :
GT-EP-0001



Description :

20MHz bandwidth, High sensitivity, Full 8 x 10cm rectangular screen with scale illumination, Excellent triggering facilities, Ideally suited for general service and education applications, A well presented front panel ensures quick and easy operation, Variable hold-off function, Convenient X-Y operation mode allows phase difference between measurements between two waveforms, Carrying handle, which may be locked into any position, also doubles as a tilt stand, Supplied complete with mains lead and comprehensive instruction manual.

NAUGRA[®]

School Educational Instruments

Website: [www.https://www.schooleducationalinstruments.com/](https://www.schooleducationalinstruments.com/), **Email:** sales@schooleducationalinstruments.com

Address: Ambala Cantt, Haryana, India **Phone:** 91-0171-2643080