

Email :
sales@schooleducationalinstruments.com
Phone: +91-0171-2601773

Product Name :
Clamp Multimeter

Product Code :
GT-EP-0004

Description :

Portable, clamp multimeter. The kit contains: Clamp Multimeter, 1x 9V PP3 battery, test leads, carry case and instructions on use. Large LCD screen 15mm x 47mm, 3½ digits LCD with a maximum reading of 1999, Highly accurate, Data hold for easy reading, Audible continuity, Low power indication

NAUGRA[®]

School Educational Instruments

Website: [www.https://www.schooleducationalinstruments.com/](https://www.schooleducationalinstruments.com/), **Email:** sales@schooleducationalinstruments.com

Address: Ambala Cantt, Haryana, India **Phone:** 91-0171-2643080